

Bodystat Product Comparison	Bodystat 1500	Bodystat 1500MDD	Quadscan 4000
Fat Mass	x	x	x
Lean / Fat Free Mass	x	x	x
Dry Lean Mass	x	x	x
Total Weight (Max 300 kg)	x	x	x
Total Body Water	x	x	x
Basal Metabolic Rate	x	x	x
BMR / Body Weight	x	x	x
Activity Energy Requirement	x	x	x
Waist:Hip Ratio	x	x	x
Body Mass Index	x	x	x
Body Fat Mass Index		x	x
Fat Free Mass Index		x	x
WELLNESS MARKER		x	
Extra-cellular Water			x
Intra-cellular Water			x
Nutrition Index			x
Body Cell Mass			x
3rd Space Water			x
ILLNESS MARKER			
Segmental Impedance Analysis			x
Frequency Measured	50 kHz	5 / 50 kHz	5 / 50 / 100 / 200 kHz
Impedance	x	x	x
Resistance		x	x
Reactance Phase Angle		x	x
Normal Ranges	x	x	x
ADDITIONAL FEATURES			
Automatic 100 Test Memory Storage for subsequent recall	x	x	x
Batteries AA	6	6	6
Real time clock		x	x
Bluetooth Communications	Only supplied with software / printer	x	x
European Foreign Languages	x	x	x