

UniGym

ALL-IN-ONE HOME GYM SYSTEM




ONE DEVICE TO MEET ALL YOUR NEEDS

- ⊙ Lightweight & Portable, workouts anytime, anywhere
- ⊙ Quick-Switch between 6 Forms, covering upper and lower limbs joints
- ⊙ 4 Modes, 10+ Exercises, with real-time adjustable training parameters
- ⊙ Feature-Rich APP with training plans, engaging games and training reports
- ⊙ 5 Safety Measures, built-in suction cup, ensuring worry-free training

UniGym All-in-One Home Gym System

Product Introduction


UniGym All-in-One Home Gym System is the world's first device covering hand, upper, and lower limb rehabilitation. With modular design, it supports 10+ exercises for fine and gross motor skills, featuring 3 training modules and 5 accessories. Lightweight and portable, its electric suction cup allows secure use on any surface. Paired with a smart app, UniGym offers 4 training modes, personalized plans, real-time adjustments, video tutorials, engaging games, and detailed reports for complete home exercise.



Product Components



One Device for Full-Body Coverage



Shoulder



Elbow



Hip



Knee



Hand



Forearm



Arm



Leg

Train Anytime, Anywhere



Lightweight and Portable Flexible Training

Weight: 6 kg Bottom Diameter: 22 cm
Battery Life: 12 hours




Tabletop



Door



Floor




Wall

Basic Edition



Includes: Main Unit, Rotation Module, App

Full Edition



Includes: Main Unit, Rotation Module, Pull Rope Module, Pneumatic Glove, App


Kids Edition




Includes: Main Unit, Rotation Module, Pull Rope Module, Pneumatic Glove, App

Smarter and More Fun

User-Friendly Interface




Training Parameters



Training Plan



Training Data



3D Racing



Baseball Battle



Home Defense



Travel Projector



More Stable, More Secure




Stable and Secure During Training

Suction Cup Resists
20 kg Lateral Pull
45 kg Vertical Pull



5


Safety Protection




Spasm Detection



Torque Limit



Speed Limit



Device Pose Detection



Emergency Stop

Specification

Dimensions	H10.7in*W19.4in*L8.7in (H27cm*W49 cm*L22 cm)	Weight	13.2 lbs (6 kg)
Battery Capacity (Wireless)	Power Lithium Battery 6400mAh	Battery Life (Wireless)	12h
Charging Voltage (Wireless)	110~220 V AC	Working Voltage (Wireless)	110~220VAC
Suction Cup Diameter	5.8 in(14.7 cm)	Suction Cup Holding Force	Lateral 44.1lbs (20kg) Longitudinal 99.2lbs(45kg)
Pull Rope Length	0-55.1in(0-140 cm)	Pull Rope Resistance	2-10kg
Rotation Diameter	7 in(17.8 cm)	Rotational Torque	0-6 Nm
Rotational Speed	0-60 rpm	APP OS	iOS/Android

May the Health and Happiness be with Everyone
